


How to measure for Venoflex Kokoon Knee High Socks

The best time to measure is first thing in the morning

How to measure for below knee socks

1. **Ankle (Cb)** measure the circumference at the ankle just above the bone
2. **Calf (Cc)** measure the circumference at the widest part of the calf

Length (CD) measure from the heel at the ground to just below the knee cap



Size	cB (cm)	cC (cm)
1	19-21	30-35
2	21-24	32-38
3	24-26	34-42
4	26-29	37-46
5 *	29-31	43-49

Length	Normal	Long
	37-40cm	greater than 40cm

Available in Compression 1 (10-15mmHg) and Compression 2 (20-36mmHg)

*** Size 5 is only available in Compression 1**